

# For Research Use Only (RUO)

# Alpha-1-Antitrypsin (AAT)

Rabbit Polyclonal Antibody

# [Catalog Number]

**REF** 0012

## 【Package Size】

Ready to use:	□1mL	□2mL	□3mL	. □5mL	□6mL
Concentrated:	□0.1m	nL □0.2	mL 🗆	0.5mL [	□1.0ml
_	_				

### [Intended Use]

Rabbit Polyclonal anti-Alpha-1-Antitrypsin (AAT) antibody is intended for use to qualitatively identify AAT antigen by light microscopy in sections of formalin-fixed, paraffin-embedded tissue using IHC detection methodology.

This antibody is for research use only (RUO).

# **Specimen Collection and Preparation for Analysis**

Formalin-fixed, paraffin-embedded tissues.

Each section should be cut to the appropriate thickness  $(2-5 \mu m)$  for the primary antibody being used and placed on a positively charged glass microscope slide.

#### **Storage and Handling**

Store at 2-8°C. Do not freeze.

Do not use product beyond the expiration date for the storage method.

#### Reagents Provided

Buffer: 10mM pH 7.4 Phosphate Puffer Saline (PBS).

Stabilizer: 0.05% bovine serum (BSA).

Preservative: 0.05% sodium azide (NaN<sub>3</sub>).

Ready-to-use antibody concentration: 2-5µg/mL.

Concentrated antibody concentration: 50-200µg/mL.

### **Staining Procedure**

 Deparaffinized slides in 3 changes of xylene (or Dewax solution), 10 minutes each. and hydrate through a graded series alcohols.

- 2. Wash the section in 90%, 80% and 70% ethyl alcohol for 10 minutes each.
- 3. Rinse in distilled water, 2 x 5 minutes.
- 4. Block the endogenous peroxidase by incubating the tissue in 3% hydrogen peroxide  $(H_2O_2)$  for 10 minutes.
- 5. Wash in distilled water, 2 x 5 minutes.
- Antigen retrieval: Place slides in a pressure cooker filled with Epitope Retrieval Solution (Citrate, pH 6.0) buffer.
- 7. Wash in PBS 2 x 5 minutes.
- 8. Concentrated Antibody Dilution

Suggested Dilution: 1:50-1:100

The optimal dilution for a specific application under a given set of experimental conditions should be determined by the investigator.

- 9. Add 100 $\mu$ L primary antibody, Incubate for 30 minutes . Wash in PBS 2 x 5 minutes.
- 10. Add 100μL secondary antibody (the protocol depends on the supplier), and proceed to standard immunohistochemistry protocol. Wash in PBS 2 x 5 minutes.
- Add 100μL DAB solution (the protocol depends on the supplier), Incubate for 2-10 minutes. Wash in PBS 2 x 5 minutes.
- Counterstain with hematoxylin. Rinse with deionized water.

# [Contact Information]



## Shanghai Long Island Antibody Diagnostica Inc.

Add: 100A, Building #1, No. 288 Hangnanzhi Road,

Zhuanghang Industrial Zone, Fengxian District, 201415

Shanghai, People's Republic of China

Tel: +86 21-64910505; 400-920-0015

Fax: +86 21-57469228

Email: sales@longislandab.com

Web: www.longislandab.com

Revision: 1.0 1/1 Date: 2020.07.01