

For Research Use Only (RUO)

Smooth Muscle Myosin Heavy Chain (SM-MHC)

Mouse Monoclonal Antibody

【Catalog Number】

REF	0413
-----	------

【Package Size】

Ready to use: \Box 1mL \Box 2mL \Box 3mL \Box 5mL \Box 6mL
Concentrated: \square 0.1mL \square 0.2 mL \square 0.5mL \square 1.0mL
7 1.4.1.1.1. 7

[Intended Use]

Mouse Monoclonal anti-Smooth Muscle Myosin Heavy Chain (SM-MHC) antibody is intended for use to identify SM-MHC antigen by light qualitatively formalin-fixed. microscopy in sections paraffin-embedded tissue usina IHC detection methodology.

This antibody is for research use only (RUO).

[Specimen Collection and Preparation for Analysis]

Formalin-fixed, paraffin-embedded tissues.

Each section should be cut to the appropriate thickness (2-5 μ m) for the primary antibody being used and placed on a positively charged glass microscope slide.

Storage and Handling

Store at 2-8°C. Do not freeze.

Do not use product beyond the expiration date for the storage method.

[Reagents Provided]

Clone: MYH11/923

Buffer: 10mM pH 7.4 Phosphate Puffer Saline (PBS).

Stabilizer: 0.05% bovine serum (BSA). Preservative: 0.05% sodium azide (NaN₃).

Ready-to-use antibody concentration: 2-5µg/mL.

Concentrated antibody concentration: 50-200µg/mL.

Staining Procedure

 Deparaffinized slides in 3 changes of xylene (or Dewax solution), 10 minutes each. and hydrate through a graded series alcohols.

- 2. Wash the section in 90%, 80% and 70% ethyl alcohol for 10 minutes each.
- 3. Rinse in distilled water. 2 x 5 minutes.
- 4. Block the endogenous peroxidase by incubating the tissue in 3% hydrogen peroxide (H_2O_2) for 10 minutes.
- 5. Wash in distilled water, 2 x 5 minutes.
- Antigen retrieval: Place slides in a pressure cooker filled with Epitope Retrieval Solution (Citrate, pH 6.0) buffer.
- 7. Wash in PBS 2 x 5 minutes.
- 8. Concentrated Antibody Dilution

Suggested Dilution: 1:50-1:100

The optimal dilution for a specific application under a given set of experimental conditions should be determined by the investigator.

- 9. Add 100 μ L primary antibody, Incubate for 30 minutes . Wash in PBS 2 x 5 minutes.
- 10. Add 100μL secondary antibody (the protocol depends on the supplier), and proceed to standard immunohistochemistry protocol. Wash in PBS 2 x 5 minutes.
- Add 100μL DAB solution (the protocol depends on the supplier), Incubate for 2-10 minutes. Wash in PBS 2 x 5 minutes.
- 12. Counterstain with hematoxylin. Rinse with deionized water.

[Contact Information]



Shanghai Long Island Antibody Diagnostica Inc.

Add: 100A, Building #1, No. 288 Hangnanzhi Road,

Zhuanghang Industrial Zone, Fengxian District, 201415

Shanghai, People's Republic of China Tel: +86 21-64910505: 400-920-0015

Fax: +86 21-57469228

Email: sales@longislandab.com

Revision: 1.0 1/1 Date: 2020.07.01



Web: www.longislandab.com